

Samsung Medison is a manufacturer of medical devices. Founded in 1985, the company sells cutting-edge medical devices, including diagnostic ultrasound, digital X-ray, and blood analysis machines, and is known around the world for its R&D capabilities and advanced technologies. In 2011, it became an affiliate of Samsung Electronics.

CT-RS85 V1.0-Eda-171110-EN

- * This product, along with its various features, options and transducers, is not currently available in all countries.
- Due to regulatory reasons, its future availability cannot be guaranteed. Please contact your local sales network for further details.
- * S-Vision™ is the name of Samsung's ultrasound imaging technology.
- * S-Vue™ is the name of Samsung's advanced transducer technology.
- * S-Detect™ for Breast and S-Detect™ for Thyroids are not available in Canada.
- * Recommendations about whether results are benign or malignant are not applicable in the United States.
- * Strain value for ElastoScan™ are not applicable in Canada or the United States.
- * All clinical images in this catalog were captured using the RS80A V4.00 ultrasound system. (except S-Shearwave Imaging™_liver image)

Empowering Professionals

Ultrasound System RS85



Scan code or visit
www.samsunghealthcare.com
to learn more



SAMSUNG MEDISON CO., LTD.

© 2017 Samsung Medison All Rights Reserved.
Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.

EXPERIENCE
A New Healthcare
Solution

SAMSUNG

A New and Outstanding Experience in Diagnosis

Beyond Experience, an integrated solution engineered to offer medical professionals a new and outstanding experience in diagnosis, delivers enriched views, advanced intelligence, streamlined workflow, and patient-centered care.

RS85 is Samsung's premium ultrasound system that adopted the integrated solution. Built with exquisite image quality and expert tools, it empowers professionals to make faster and more confident decisions.



Ultrasound System **RS85**

BEYOND EXPERIENCE™

Samsung's commitment to supporting confident decision making

Advanced Intelligence



Enriched View



BEYOND EXPERIENCE™




Streamlined Workflow

Patient-centered Care




Enriched View Exquisite Image Quality




S-Vision™ Imaging Engine

With the S-Vision™ imaging engine built into RS85, the digital signals produce clear, detailed resolution and tissue uniformity for various types of applications in general imaging.



S-Harmonic™

This new harmonic technology improves image clarity, near to far. Reducing signal noise, S-Harmonic™ provides more uniform ultrasound images. Combined with the S-Vue™ transducers, S-Harmonic™ takes RS85 image quality one step further.




Gallbladder

Improved Diagnostic Confidence

Samsung's image enhancing and artifact suppressing technologies and S-Vue™ transducers together provide clear, detailed imaging that you can count on to help improve diagnostic confidence and imaging continuity.


S-Vue™ Transducers

S-Vue™ transducers provide more efficient piezoelectric properties, resulting in wider bandwidths that enable better penetration and higher quality resolution on even challenging patients.



HQ Vision™


HQ Vision™ is a new, advanced technology for visualizing anatomical structures. With improved image clarity, this feature helps make a reliable diagnosis quickly.



Knee

ClearVision

The noise reduction filter improves edge enhancement and creates sharper 2D images for optimal diagnostic performance. The integration of specialized Samsung technology results in a notable improvement in image quality. In addition, ClearVision provides application-specific optimization and advanced temporal resolution in live scan mode.



Abdomen




Enriched View
Expert Tools

More Valuable Information


Expert tools offer new perspectives and provide additional information for confident decision making.

MV-Flow™

MV-Flow™ offers a novel alternative to power Doppler for visualizing slow flow microvascularized structures. High frame rates and advanced filtering enable MV-Flow™ to provide a detailed view of blood flow in relation to surrounding tissue or pathology with enhanced spatial resolution and temporal resolution.



Kidney



Thyroid

CEUS+

CEUS+ technology uses the unique properties of ultrasound contrast agents. When stimulated with low acoustic pressure, the oscillating microbubbles reflect both fundamental and harmonic frequency signals. In addition, Samsung's technologies, VesselMax™ and FlowMax™, provide a clear visualization of vessels and blood flow for a more informed and confident diagnosis.

S-Fusion™

S-Fusion™ enables simultaneous localization of a lesion using real-time ultrasound in conjunction with other volumetric imaging modalities. Samsung's Auto Registration helps quickly and precisely fuses the images, increasing efficiency and reducing procedure time. S-Fusion™ enables precise targeting during interventional and other advanced clinical procedures.

S-Fusion™ for Prostate

S-Fusion™ for Prostate allows precise targeting during prostate biopsies. Based on 3D models created with MR data sets, S-Fusion™ for Prostate provides biopsy guidance to help safely navigate and target the prostate.




Advanced Intelligence



Arterial Analysis™


Arterial Analysis™ detects functional changes of vessels, providing measurement values such as the stiffness, intima-media thickness and pulse wave velocity of the common carotid artery. Since the functional changes occur before morphological changes, this technology supports the early detection of cardiovascular disease.



Arterial Analysis™

S-3D Arterial Analysis™

S-3D Arterial Analysis™ simplifies volume measurement of arterial plaque, providing 3D vessel modeling. With Samsung's S-3D Arterial Analysis™, obtaining information on the arterial plaque volume is surprisingly fast and easy even on difficult patients. In addition, it allows you to track the morphological changes of the artery.




S-3D Arterial Analysis™

Increased Consistency

With its advanced intelligent solutions, including an extensive range of quantification functions, RS85 provides measurement consistency while reducing variability between users.

S-Shearwave Imaging™

S-Shearwave Imaging™ allows for non-invasive assessment of the stiffness of tissue/lesions in the breast and liver, by providing an advanced level of diagnostic information. The color-coded elastogram, quantitative measurements (in kPa or m/s), dual or single display option, and user-selectable ROI (position and size) functions are especially useful for the accurate diagnosis of breast and liver diseases.




Liver

Breast

S-Detect™ for Breast

S-Detect™ for Breast helps standardize reporting and classification of suspicious breast lesions by incorporating BIRADS @ ATLAS* (Breast Imaging-Reporting and Data System, Atlas) into the tool. When the user selects a region of interest, S-Detect™ for Breast automatically defines the lesion boundaries, provides lexicon classification options, and images export for an enhanced and streamlined workflow.

*Registered trademark of the American College of Radiology. All rights reserved.




Breast

S-Detect™ for Thyroid

S-Detect™ for Thyroid uses the advanced technology based on K-TIRADS, RUSS and ATA guideline* in detecting and classifying suspicious thyroid lesions semi-automatically. This technology helps you diagnose your patients with confidence and ease, providing accurate, consistent results and an automatic reporting feature.

* K-TIRADS: Korean-Thyroid Imaging Reporting and Data System
RUSS: Russ' TIRADS
ATA: American Thyroid Association




Thyroid


Enhanced Efficiency

The RS85 has been designed to streamline your workflow by enhancing efficiency through reducing keystrokes and by combining multiple actions into one.

Quick Preset
With one touch, the user can select the most common transducer and preset combinations. Quick Preset maximizes efficiency to make a full day of scanning simple and easy.



Touch Customization
Samsung has made a customizable touchscreen interface that allows the user to move frequently used functions to the first page, keeping the focus on the patient instead of the system.



13.3-inch Tilting Touch Screen
Samsung's tilting touch screen can be adjusted to accommodate any user's viewing preferences within any scanning environment.



6-way Control Panel
The RS85's 6-way adjustable control panel optimizes your work environment to reduce repetitive motions stress. When it's in off-mode, the control panel returns to the home position, allowing for easier and enhanced mobility.



Comfort and Safety


The patient-centric gel warmer and low noise functions work together to help your patient feel more comfortable while they're being examined.



Low Noise Levels
The RS85 is a quiet device that allows physical exams to be performed, including auscultation, while the ultrasound system is turned on.



Gel Warmer
Samsung's two-level adjustable gel warmer keeps ultrasound gel at a comfortable temperature.



Central Lock
A single pedal controls a central lock mechanism to conveniently secure the console in place. This results in more efficient movements while the user is performing scanning procedures.




Maneuverable Wheel
4 swivel wheels allow easy steering, and a locking function.



Image Gallery




Bowel




Liver



Pancreas




CCA




E-Strain™




Median nerve




Auto IMT+




Kidney




Breast mass



E-Thyroid™



Strain+



4 Chamber

Comprehensive Selection of Transducers

S-Vue™ transducers

Curved array transducers



Application: abdomen, obstetrics, gynecology



Application: abdomen, obstetrics, gynecology, contrast



Application: abdomen, obstetrics, gynecology

Volume transducer

Linear array transducers



Application: small parts, vascular, musculoskeletal



Application: small parts, vascular, musculoskeletal



Application: small parts, vascular, musculoskeletal



Application: small parts, vascular, musculoskeletal



Application: small parts, vascular, musculoskeletal, abdomen



Application: musculoskeletal

Curved array transducers



Application: abdomen, obstetrics, gynecology



Application: pediatric, vascular, gynecology

Volume transducers



Application: musculoskeletal, small parts, vascular



Application: obstetrics, gynecology, urology

Endo-cavity transducer



Application: obstetrics, gynecology, urology

Phased array transducers



Application: cardiac, TCD, abdomen



Application: cardiac, pediatric, abdomen



Application: cardiac, pediatric

CW transducers



Application: cardiac



Application: cardiac



Application: cardiac

TEE transducer



Application: cardiac

* Some of the transducers may not be available in some countries.